



towel & tissue

THE J2 WAY



Major Cost Factors

Labor: In a busy facility, the labor cost of refilling towel and tissue dispensers could easily cost thousands of dollars per year - the single largest expense for towel and tissue.

Usage: Excessive usage is the second highest cost factor for towel and tissue.

Waste: Towel and tissue can be wasted due to stub roll waste, moisture, theft, and/or vandalism to product or dispensing systems.

Product: Some people will try to save money by buying lighter weight products such as 1-ply tissue paper, only to discover that people will use more of it - negating any savings. Others will shop for the lowest “case cost” only to discover that it is actually more expensive by the foot or by the use.

Batteries: Electronic dispensers are excellent in certain locations such as healthcare or food service but are not necessary nor cost-effective in most other environments due to battery costs which can cost hundreds or thousands per year.



The lack of usable paper is the most common restroom complaint

which explains why maintenance personnel will discard 20% of useful product to reduce the chance of run-outs.

WHAT IS THE J2 WAY?



- Reduce customer costs with better solutions
- Improve customers' productivity
- Make acquiring products easier for our customers
- Protect our customers' brand

“We promise to do our best each and every time to meet our customers' needs with exceptional service and solutions.”



Mechanical Roll
Towel Dispenser
T880B



Mini Smart Core Twin
Tissue Dispenser
T200B



Mini Smart Core Twin
Tissue Dispenser
T400B



3 Roll Carousel Smart
Core Tissue Dispenser
3342



Jumbo Smart Core Twin
Tissue Dispenser
T370B

By switching from standard bath tissue to a Transcend® Twin Tissue System, labor is reduced by an amazing 89.8%!

Solutions

High Capacity: Greatly reducing labor costs is the main benefit of high-capacity dispensing systems.

Auto-Cut Dispensing Systems: Auto-cut dispensers regulate and help control usage by automatically cutting the towel at a predetermined length.

Locking Mechanisms: A secure locking mechanism reduces costs (and maintenance) by ensuring that the paper products stay where they belong and are protected from damage due to vandalism.

Secure Dispensing Systems: Dispensers are manufactured using impact-resistant plastic to help withstand abuse and vandalism, and drainage channels are designed to move moisture away from the paper products inside the dispensers, resulting in less waste.

Mechanical Dispensing Systems: Dispensers are manufactured using impact resistant plastic to help withstand abuse and vandalism, and drainage channels are designed to move moisture away from the paper products inside the dispensers, resulting in less waste.

Quality Assurance: Tissue and towels made to a quality standard that ensures the quality and quantity of products are consistent and accurate.

Replacement: J2 will fix or replace damaged dispenser free of charge.



TOWEL & TISSUE PRODUCTS



| NAME | | ITEM NO. | UNIT | DESCRIPTION | PRICE |
|--|---|----------|------|--|-------------------------------------|
| Mechanical Roll Towel Dispenser |  | T880B | EACH | Max. Capacity: 1,000' Avg. Hand-Dries: 545 ADA Compliant: Yes | Use: 813BP or 813NP Towels |
| Electronic Roll Towel System |  | 8800EP | EACH | Max. Capacity: 1,000' Avg. Hand-Dries: 545 ADA Compliant: Yes | Use: 813BP or 813NP Towels |
| White Roll Towels |  | 813BP | CASE | 1,000 feet per roll 6 rolls per case Green Seal™ Certified | |
| Natural Roll Towels |  | 813NP | CASE | 1,000 feet per roll 6 rolls per case Green Seal™ Certified | |
| Mini Smart Core Twin Tissue Dispenser |  | T200B | EACH | Max. Capacity: 840' Avg. Uses: 260.5 ADA Compliant: Yes | Use: PR420 or PR420VT Tissue |
| Transcend 2-ply Tissue |  | PR420 | EACH | 420 feet per roll 24 rolls per case Green Seal™ Certified | |
| Jumbo Smart Core Twin Tissue Dispenser |  | T370B | CASE | Max. Capacity: 2,290' Avg. Uses: 710 ADA Compliant: Yes | Use: 1145 Tissue |
| Transcend 2-ply Jumbo Tissue |  | 1145 | CASE | 1145 feet per roll 12 rolls per case Green Seal™ Certified | |
| Transcend 2-ply Premium Jumbo Tissue |  | 1145VT | CASE | 1145 feet per roll 12 rolls per case Green Seal™ Certified | |



World Health Organization

In a comparison of methods to test efficiency of hand drying for removing bacteria from washed hands, warm air drying performed worse than drying with paper towels. Furthermore, air dryers may be less recommendable because of longer time needed to achieve dry hands, with a possible negative impact on hand hygiene compliance, and because of the aerosolization of waterborne pathogens. Ideally, drying of hands should be done by using individual paper towels.

J2 Business Products

KALISPELL
406.752.8520

support@j2op.com

shopj2.com

SALES REP:

PHONE: